


ANNUAL QA

PRECISE, EFFICIENT AND
BUILT FOR THE FUTURE



DOSEVIEW™ 3D

DOSEVIEW™ 3D

COMMITTED TO INNOVATION AND QUALITY PATIENT OUTCOMES

Making sure your patients have great treatments begins with a finely-tuned linac, and that starts with the DoseView 3D.

It isn't just another water phantom.



BEN ROBISON

MEDICAL PHYSICIST

PROVISION RADIATION THERAPY

"The DoseView 3D is an outstanding scanning tank. Setup is simple and easy. Hopefully you are only using your scanning tank once a year. DV3D has made the set up so simple that after 12 months of not using the tank it still feels familiar."

BETTER HARDWARE

ROBUST CONSTRUCTION PROVIDES THE FOUNDATION FOR BETTER DATA

STRONG

Manufactured with a rigid, aluminum frame and stainless steel leadscrews, the DoseView 3D provides a consistent platform for repeated QA testing. Meticulous construction ensures consistent measurement accuracy within ± 0.1 mm per axis.

ACCURATE

The Cartesian coordinate system and ability to set isocenter to any point within the tank avoids rotational errors possible with cylindrical tanks.



CONVENIENT

The wireless pendant with audible and visual cues and a backlit display make configuring and controlling the DoseView 3D easy to perform and remember.

EASY

Three leveling screws are mounted above the phantom's scanning volume, making it easy to fine-tune even when the phantom is filled with water.

AUTOMATIC
DETECTOR ALIGNMENT
SYSTEM

Swap detectors without resetting isocenter/origin with the Detector Alignment System. The system aligns the center of each detector.



SUPPORTS
VARIAN® HALCYON™

BUILT FOR ACCURACY




HIGH-SPEED ELECTROMETER

The DoseView 3D Electrometer is a dual-channel, fast-acquisition electrometer that continues Standard Imaging's renowned reputation for precision and efficiency. Bias voltage can be enabled independently per channel for chamber/diode simultaneous operation. Low noise performance helps ensure post-processing operations are kept to a minimum.

ACCURATE, RELIABLE EXRADIN ION CHAMBERS


For over 40 years, Exradin's quality has been recognized by the top research institutes and standards laboratories around the world. The DoseView 3D comes with two Exradin A28 Ion Chambers (0.125cc), producing sustained accuracy and confidence in your QA testing.



VARIABLE STEP SIZE ACQUISITION

The DoseView 3D can be programmed to move in varying step sizes throughout a scan. Implement tighter steps toward the penumbra to account for increased variability without significantly increasing overall scan time.

TIME SAVING FEATURES



EXPEDITE ACCURATE SET UP

The Automatic Detector Alignment System places the center of all detectors at the same position relative to the water's surface, allowing users to swap detectors without resetting isocenter/origins. Initial setup is performed with the crosshair alignment jig. This bracket system is compatible with Exradin Ion Chambers and supports chambers and diodes from most other manufacturers. This system ensures faster setup time and provides confidence in accurate leveling and beam center positioning.

AUTOMATICALLY FIND BEAM CENTER

Once leveling is performed, the DoseView 3D determines the center of a radiation field by locating the field edges using an ion chamber or diode. This test can be performed at multiple depths to identify any issues with system leveling and/or gantry alignment.

INCREDIBLE, ALL NEW SOFTWARE DESIGN

THE EASIEST, MOST INTUITIVE WATER PHANTOM SOFTWARE AVAILABLE

- A wireless connection to the hardware provides easy & convenient control
- Import scans from your previous water scanning systems for comparison & continuity.
- Accuracy of system setup is confirmed with automatic checks that ensure detector centering and gantry angle positioning.
- Scan queue setup and measurement is automated eliminating repetitive tasks and saving time.
- Scan speed and step size automatically optimized according to dose rate.
- Guided beam scanning performing next sequential group of measurements in queue.
- SmartSelect detects scans outside queue parameters and automatically deselects them.
- Simple, yet powerful scan queue creation and editing.
- Easy export to TPS and CSV

QUEUES

A queue creation wizard helps you choose your scanning modality (photons, electrons or flattening filter free photons), modifiers, energies, field sizes and scan types. All combinations of your choices are created. No copying and editing individual scan groups. Make your scan queues in a matter of minutes...or less. Continuous scanning and multiple scan zones are available at no extra cost.

Add Queue

General

6x Annual Scans, Scans, Water, Jaws

Complete

Energies

6 MV

Complete

Field Sizes

20 x 20, 10 x 10, 5 x 5

Complete

Scan Types

CrossPlane, InPlane, Depth

Complete

Depths

Not Set

CAX Offsets

Not Set

Scan Mode

Complete

SETUP

The Measure page has a setup wizard that takes you through the process of getting your tank ready to use. When all expander title bars have results in blue, you are ready to scan.

Hardware

Queue

Machine

SI Irradiator, 60, 80

Dose Rates

Configured

Water Phantom

Ready

DEMO

Connect

Initialize

Detectors

Exradin A28, Exradin A4

Field:

Exradin A28

Reference:

Exradin A4

Do not use a reference detector

Electrometer

Field: --, Ref: --

Position

CP: 0.00 cm, IP: 0.00 cm, Depth: 0.00 cm


Phantom Orientation

Controller Right




MEASURE

Scan queues are loaded into a queue runner. As scan data is acquired, it is immediately stored in a database for future access, analysis, and export. Raw scan data and the methods used to modify it are stored in the database, so you always know what your raw data looked like and always know how you processed it for your needs.



ANALYZE

Comprehensive scan processing and analysis tools include profile centering, mirroring, symmetric averaging, shifting & normalization of profiles & depth doses, smoothing with boxcar & gaussian filters, and dose conversion using TG-51 and TRS-398.



LIFT AND RESERVOIR CARTS

FULLY INTEGRATED WITH LIFT AND RESERVOIR IN ONE

The DoseView 3D's cart contains both an electronic lift mechanism and 60 gallon (265 liter) water reservoir, resulting in hassle-free storage and setup. Additionally, a convenient storage area holds the power supply and other accessories. The lift and pump are powered by a single power cable, and a power pass-through cable allows convenient connection to the water phantom.



EXTENSIVE ADJUSTMENTS CAPABILITIES

The Precision Positioning Platform provides an ideal medium for maneuvering the DoseView 3D. This added flexibility allows subtle X and Y axis movements up to ±12.5mm, phantom rotation of ±1° and positive engagement at 10°, 45°, and 90° intervals. Achieve consistent detector orientation during in-plane, cross plane and diagonal scans.

EFFICIENCY LIFT CART

An alternate, portable lift cart is also available. The electric lift table has a capacity of 900 pounds and is capable of over 15 inches of vertical travel.



24/7 SUPPORT

Standard Imaging's acclaimed customer service and support ensures an easy transition of the DoseView 3D into your QA workflow and will guide you through any issues that may arise in the future. Standard Imaging also provides 24/7 phone support for DoseView 3D in USA.

For more information on warranties and additional maintenance and support packages please contact your Standard Imaging Account Manager or regional distributor.

DOSEVIEW 3D SPECIFICATIONS

MOTION CONTROL SYSTEM	
TRAVEL SPEED — 50 mm/s // POSITIONING ACCURACY — ± 0.1 mm per axis // POSITIONING REPEATABILITY — ± 0.1 mm per axis	
PC COMMUNICATION — Wireless or wired via RS-232 // CONTROL METHOD — PC or via wireless pendant	
WATER PHANTOM (ACRYLIC TANK ONLY)	
OUTER DIMENSIONS — [length × width × height] 683 mm x 692 mm x 542 mm	
SCANNING DIMENSIONS — [length × width × height] 480 mm x 480 mm x 410 mm	
WALL THICKNESS — 19 mm // Replaceable Fill/Drain Port	
LIFT CART AND RESERVOIR	
OUTER DIMENSIONS — [length × width] 1247 mm x 762 mm // VERTICAL RANGE — 685 mm – 1185 mm (tank base to floor)	
WATER PUMP — Electric fill, gravity drain // WATER CAPACITY — 60 gal (227 liters)	
FILL SPEED — 6-7 min // DRAIN SPEED — 16-20 min	
PRECISION POSITIONING PLATFORM	
X / Y FINE ADJUSTMENT — ±12.5 mm // FINE ROTATIONAL ADJUSTMENT — ±1°	
DISCREET ENGAGEMENT — 10°, 45°, and 90° intervals	
DOSEVIEW 3D ELECTROMETER	
CHANNELS — 2 // BIAS VOLTAGE — 0, ±150 to 450 (VDC) in 50 volt increments // RANGE — 2 pC – 999,999 nC	
RESOLUTION — 10 fC // CONNECTOR TYPE — Triaxial BNC or TNC (BNC unless specified)	
SOFTWARE/COMPUTER REQUIREMENTS	
OPERATING SYSTEM — Windows® 10 Professional, 64-bit Recommended	
RUNTIME ENVIRONMENT — .NET 4.5.2 // PROCESSOR — Dual Core, 1 GHz; Quad Core, 2 GHz Recommended	
MEMORY — 32-bit OS: 2 GB, 4 GB Recommended 64-bit OS: 4 GB, 8 GB Recommended	
HARD DRIVE — 1 GB free space for initial software installation. Typically 2.5 MB/Year/Linac disk space growth. 25% free space Recommended.	
SCREEN RESOLUTION — 1440 x 900 or greater // CD-ROM DRIVE — 2X speed or greater	
CONNECTIVITY — 9-pin serial RS-232 port or USB port with USB to RS-232 adapter	

Microsoft® and Windows® are registered trademarks of Microsoft Corporation. Pinnacle3® is a registered trademark of Philips Electronics N.V. Eclipse™ and Halcyon™ are trademarks of Varian Medical Systems Inc. Monaco® is a registered trademark of Elekta AB.

Specifications subject to change without notice.

DoseView 3D REF 92260, DoseView 3D Lift and Reservoir Cart REF 72260

ANNUAL QA

STANDARD IMAGING®



www.standardimaging.com

800-261-4446 . PH 608-831-0025 . FAX 608-831-2202

3120 Deming Way Middleton WI 53562 USA

© 2020 Standard Imaging, Inc.

1322-29 (9.20)

